

**ANSI-Dsen-Thermoelement, glasfaserisoliertes
Kabel mit Edelstahl-bergeflecht, das in blanken
Schwunzen endet, Miniatur- oder Standardstecker -
Typ K**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

General Description	Thermocouple housed in a stainless steel threaded bolt
Thermocouple Type	K
Thermocouple Junction	Grounded hot junction for fast response, M6x1mm thread
Cable Length	1, 2 or 5m
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Bare tails, Miniature or Standard Plug (with cable clamps)
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	ANSI MC96.1