



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI-Unterlegscheiben-Thermoelement, glasfaserisoliertes Kabel mit Edelstahl-bergeflecht, das in blanken Schwunzen oder Standardstecker endet Typ K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

General Description	Thermocouple housed in a washer eyelet
Thermocouple Type	K
Thermocouple Junction	Grounded tip for quick response
Cable Length	2m
Cable Type	Glassfibre insulated, stainless steel overbraided
Core / Strands	7/0.2mm
Termination Type	Either tails or standard plug termination with cable clamps
Max. Temperature	+350°C
Standards Met	ANSI
Washer	3.5mm, 4.0mm, 5.0mm, 6.0mm, 8.0mm & 10.0mm clearance hole available