



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Unterlegscheiben-Thermoelement, glasfaserisoliertes Kabel mit Edelstahl-bergeflecht - Typ K JIS mit blanken Endverschlüssen und 10,0 mm Unterlegscheibenlochgre



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0981-001
General Description	Thermocouple housed in a washer eyelet
Thermocouple Type	K JIS
Thermocouple Junction	Grounded tip for quick response
Cable Length	2m
Cable Type	Glassfibre insulated, stainless steel overbraided
Core / Strands	7/02.mm
Termination Type	Bare Tails
Max. Temperature	+350°C
Standards Met	JIS (Japanese Industrial Standards)
Washer	10.0mm clearance hole