


# Standard-Thermoelement-Steckverbinder fr die Schalttafelmontage mit Edelstahlhalterung IS-K-SSPF Typ K IEC



**Standard in-line Socket with stainless steel mounting bracket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Standard-Thermoelement-Steckverbinder fr die Schalttafelmontage mit Edelstahlhalterung IS-K-SSPF Typ K IEC**

Vorgesehen fr den Einsatz in eher industriellen Prozessen, kann aber beliebig in allen Anwendungen eingesetzt werden

**Leistungsbeschreibung****Specifications**

<b>Product Code</b>	XE-1408-001
<b>General Description</b>	Intended for use in more industrial processes but can be used in all applications as desired
<b>Thermocouple Type</b>	K IEC
<b>Connector Type</b>	Panel Mount with SS Bracket
<b>Connector Size</b>	Standard
<b>Colour</b>	Green
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS