



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

PTFE-Oliven



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

PTFE-Oliven. Erhältlich für Sondendurchmesser von 1,5 mm bis 8,0 mm

Specifications

General Description	Olives
Material	PTFE
Probe Diameter	Suitable for probe diameters between 1.5mm & 8.0mm