


Standard Panel Mount Thermoelement Connector mit Edelstahlhalterung IS-E-SSPF Typ E IEC



Standard in-line Socket with stainless steel mounting bracket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard Panel Mount Thermoelement Connector mit Edelstahlhalterung IS-E-SSPF Typ E IEC

Für den Einsatz in industriellen Prozessen vorgesehen, kann aber in allen Anwendungen nach Belieben eingesetzt werden

Spezifikationen**Specifications**

| | |
|----------------------------|--|
| Product Code | XE-1638-001 |
| General Description | Intended for use in more industrial processes but can be used in all applications as desired |
| Thermocouple Type | E IEC |
| Connector Type | Socket with SS Bracket |
| Connector Size | Standard |
| Colour | Violet |
| Max. Temperature | 220°C |
| +Positive Contact | Nickel Chromium |
| -Negative Contact | Copper Nickel |
| Standards Met | BS EN 50212:1997. RoHS |