


PFA-isoliertes, doppelt verdrilltes IEC-Thermoelement mit freiliegender Sperrschicht und Bare-Tails - Typen K, J, T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermoelement Typ K, J oder T. Schnell ansprechende, geschweifte, freiliegende Verbindung in Kabelspezifikation der Klasse 1 gem BS EN 60584-3:2008 mit Bare-Tails-Anschluss.

Ideal für Test- und Entwicklungsanwendungen

Leistungsbeschreibung

Specifications

General Description	Fast response, welded exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	K, J & T
Cable Length	1, 2, 5 & 10m
Cable Type	Teflon® PFA insulated twin-twisted lead
Core / Strands	1/0.2mm, 1/0.315mm, 7/0.2mm, 1/0.5mm
Cable Termination	Bare Tails
Max. Temperature	+250°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC-584-3