


## **PFA-isoliertes, doppelt verdrilltes IEC-Thermoelement mit freiliegender Sperrschicht und Bare-Tails - Typen K, J, T**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermoelement Typ K, J oder T. Schnell ansprechende, geschweifte, freiliegende Verbindung in Kabelspezifikation der Klasse 1 gem BS EN 60584-3:2008 mit Bare-Tails-Anschluss.

Ideal für Test- und Entwicklungsanwendungen

### Leistungsbeschreibung

## Specifications

<b>General Description</b>	Fast response, welded exposed junction
<b>Sensor Type</b>	Exposed Junction
<b>Thermocouple Type</b>	K, J & T
<b>Cable Length</b>	1, 2, 5 & 10m
<b>Cable Type</b>	Teflon® PFA insulated twin-twisted lead
<b>Core / Strands</b>	1/0.2mm, 1/0.315mm, 7/0.2mm, 1/0.5mm
<b>Cable Termination</b>	Bare Tails
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-75°C
<b>Tolerance</b>	Class 1
<b>Standards Met</b>	IEC-584-3