

## Sonde mit flacher Spitze Typ RTD / PRT 2M Länge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Universelle Pt100-Sonde

An der Spitze abgeflacht für ein verbessertes thermisches Ansprechverhalten

2-Draht-Konfiguration

Klasse B, nach IEC 751

Mantel aus Edelstahl 316, 4,0 mm Durchmesser x 150 mm lang

1 Meter langes 7/0,2 mm Teflon-isoliertes, verdrehtes 2-adriges Kabel

Inklusive 1/8" BSPP Messing-Klemmverschraubung mit Kontermutter

Temperaturbereich -50C bis + 250C

### Leistungsbeschreibung

## Specifications

<b>Product Code</b>	XE-3698-001
<b>General Description</b>	316 stainless steel probe, 4.0mm diameter (6mm max. across flat area of tip) x 150mm long overall
<b>Sensor Type</b>	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B
<b>Probe Temperature Range</b>	-50°C to +250°C
<b>Cable Length</b>	2M
<b>Cable Type</b>	7/0.2mm Teflon® insulated 2 core twisted leads
<b>Termination Type</b>	2-wire tails
<b>Accessories</b>	1/8" BSPP brass male compression fitting with locknut included