

Flache Sonde Typ RTD / PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Universell einsetzbare Pt100-Sonde

An der Spitze abgeflacht für ein verbessertes Wärmeverhalten

2-Draht-Konfiguration

Klasse B nach IEC 751

Mantel aus Edelstahl 316, 4,0 mm Durchmesser x 150 mm Länge

1 Meter langes, 7/0,2 mm teflonisoliertes, 2-adriges, verdrehtes Kabel

Inklusive 1/8" BSPP Messing-Klemmverschraubung mit Kontermutter

Temperaturbereich -50°C bis + 250°C

Leistungsbeschreibung

Specifications

General Description	316 stainless steel probe, 4.0mm diameter (6mm max. across flat area of tip) x 150mm long overall
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B
Probe Temperature Range	-50°C to +250°C
Cable Length	1m, 2m, 5m
Cable Type	7/0.2mm Teflon® insulated 2 core twisted leads
Termination Type	2-wire tails
Accessories	1/8" BSPP brass male compression fitting with locknut included