



Miniatur-Thermoelement-Steckverbinder mit rechteckiger Faszienbuchse BM-K-FF Typ K BS



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Rechteck-Thermoelement-Steckverbinder Blendenbuchse BM-K-FF Typ K BS

Nylon-Clip für eine schnelle und einfache Befestigung an der Platte. Alle Faszienbefestigungsbuchsen sind mittels ungleich großer Stifte polarisiert

Leistungsbeschreibung**Specifications**

Product Code	XE-1383-001
General Description	All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	K BS
Connector Type	Fascia mounting socket
Connector Size	Miniature
Colour	Yellow
Dimensions	See Schematic Drawing (Product Literature Tab)
Max. Temperature	120°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS