


Standard-Rechteckige Blenden-Thermoelement-Steckverbinderbuchse BS-K-FF Typ K BS



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Rechteck-Thermoelement-Steckverbinder Blendenbuchse BS-K-FF Typ K BS

Nylon-Clip für eine schnelle und einfache Befestigung an der Platte. Alle Fasziensbefestigungsbuchsen sind mittels ungleich großer Stifte polarisiert

Leistungsbeschreibung**Specifications**

Product Code	XE-1387-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	K BS
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Yellow
Max. Temperature	120°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS