


Standard-Thermoelement-Steckverbinder In-Line-Buchse BS-K-F Typ K BS



Standard in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Thermoelement-Steckverbinder In-Line-Buchse BS-K-F Typ K BS

Akzeptiert Standard-Rundsteckerstifte, geeignet für industrielle Anwendungen

Leistungsbeschreibung**Specifications**

Product Code	XE-1389-001
General Description	Accepts standard round plug pins, contacts are polarised to prevent incorrect connection.
Thermocouple Type	K BS
Connector Type	Socket
Connector Size	Standard
Colour	Yellow
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS