


# Barriereklemme Block Thermoelement Stecker IS-CU-BTS Typ CU



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Barriereklemme Block Thermoelement Stecker IS-CU-BTS Typ CU**

Lieferung in einzelnen, ineinandergreifenden Paaren. Anschlüsse im Thermoelementmaterial

**Spezifikationen****Specifications**

<b>Product Code</b>	XE-1944-001
<b>General Description</b>	Supplied in single, interlocking pairs. Connections in thermocouple material
<b>Thermocouple Type</b>	IEC, ANSI, JIS Type CU
<b>Connector Type</b>	Barrier Terminal Blocks
<b>Colour</b>	White
<b>Dimensions</b>	See Schematic Drawing (Product Literature Tab)
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper
<b>Standards Met</b>	BS EN 50212:1997. RoHS