



Thermoelement-Flachstecker Typ J



Thermocouple Spade Terminal

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

- Crimp-auf-Flachstecker
- Preis pro Packung mit 20 Paar
- Geeignet fr IEC-, ANSI- und JIS-Kalibriertypen.
- Hergestellt aus kompatiblen Thermoelementmaterial
- Thermoelement-Typ, der auf Klemme gestempelt ist
- Suits Labfacility & andere Sperrklemmen
- Fr Drahtgren 0,5 bis 1,0 mm

Klemmenleiste

Specifications

Product Code	XE-1985-020
General Description	pack of 20 pairs
+Positive Contact	Iron
-Negative Contact	Copper Nickel