

Miniatur-Thermoelement-Steckverbinder In-Line-Buchse AM-K-FS Typ K ANSI



Miniature in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Thermoelement-Steckverbinder In-Line-Buchse AM-K-FS Typ K ANSI

Ideal for general purpose use, as all contacts are polarised, to ensure a correct connection.
These connectors accept thermocouple cable with a diameter of up to 4 mm.

Leistungsbeschreibung**Specifications**

Product Code	XE-1132-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
Thermocouple Type	K ANSI
Connector Type	Socket
Connector Size	Miniature
Colour	Yellow
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector