

## Miniatur-Thermoelement-Steckverbinder In-Line-Buchse AM-K-FS Typ K ANSI



**Miniature in-line Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur-Thermoelement-Steckverbinder In-Line-Buchse AM-K-FS Typ K ANSI**

Ideal for general purpose use, as all contacts are polarised, to ensure a correct connection.  
These connectors accept thermocouple cable with a diameter of up to 4 mm.

**Leistungsbeschreibung****Specifications**

<b>Product Code</b>	XE-1132-001
<b>General Description</b>	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
<b>Thermocouple Type</b>	K ANSI
<b>Connector Type</b>	Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Yellow
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS