


Allzweck-Thermistorsonde mit PTFE-Leitung



[Click here to view product](#)

Product Code	XE-9840
General Description	2-wire 10K NTC thermistor housed in a stainless steel sheath
Thermistor	10K NTC Thermistor
BETA 25/85	3977K
Probe Diameter	3mm, 4mm, 4.5mm & 6mm
Cable Length	2M
Cable Type	PTFE insulated twisted leads
Core / Strands	7/0.2mm
Cable Termination	Bare Tails
Max. Temperature	+150Å°C

Min. Temperature	-50Å°C
Accuracy	+/- 0.2K (over 0-70Å°C)
Length	25mm, 50mm, 100mm & 150mm
Number of Wires	2