


Thermistor Temperatursensoren fr den Auenbereich / Khlhaus



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

2-Draht-10K-NTC-Thermistor in einem wetterfesten IP67-zertifizierten Wandgehäuse für genaue Messung der Außenlufttemperatur

Hauptmerkmale:

Widerstandswert: 10K Ohm

Typ: NTC-Thermistor (Negativer Temperaturkoeffizient), was bedeutet, dass sein Widerstand mit steigender Temperatur abnimmt.

Temperaturbereich: Kann Temperaturen von -50 C bis +150 C genau messen.

Verkabelung: 2-Draht-Konfiguration für einfache Konnektivität.

IP-Schutzart: IP67

Anträge:**HLK-Systeme**

: überwacht die Lufttemperatur, um Heizungs-, Lüftungs- und Klimaanlage effizient zu regulieren.

Leistungsbeschreibung**Specifications**

Product Code	XE-9867-001
General Description	Weatherproof IP67 rated wall mounting Thermistor for accurate outdoor air temperature sensing
Sensor Type	Single Element
Thermistor	10K NTC Thermistor
BETA 25/85	3977K
Sheath Material	Perforated 316 stainless steel probe, 12.7mm x 75mm
Core / Strands	7/0.2mm
Cable Termination	Bare Tails



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Gland	M20 cable entry gland
Max. Temperature	+150°C
Min. Temperature	-50°C
IP Rating	67
Accuracy	+/- 0.2K (over 0-70°C)
Number of Wires	2