

Lufttemperatur / Thermistorsensor für den Innenbereich



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

2-Draht-10K-NTC-Thermistor in einem Innentempersensoren für die Wandmontage zur Messung der Umgebungstemperatur (nur für den Innenbereich).

Hauptmerkmale:

Widerstandswert: 10K Ohm

Typ: NTC-Thermistor (Negativer Temperaturkoeffizient), was bedeutet, dass sein Widerstand mit steigender Temperatur abnimmt.

Temperaturbereich: Kann Temperaturen von -50 C bis +70 C genau messen.

Verkabelung: 2-Draht-Konfiguration für einfache Konnektivität.

Anträge:**HLK-Systeme**

: überwacht die Lufttemperatur, um Heizungs-, Lüftungs- und Klimaanlage effizient zu regulieren.

Leistungsbeschreibung**Specifications**

Product Code	XE-9868-001
General Description	Wall mounting indoor temperature sensor for environmental temperature measurement
Sensor Type	For indoor use only
Thermistor	10K NTC Thermistor
BETA 25/85	3977K
Housing	Moulded case L85 x W85 x H30mm
Max. Temperature	+70°C
Min. Temperature	-50°C
Accuracy	+/- 0.2K (over 0-70°C)



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Features #1	Vented construction
Features #2	Surface or standard electrical conduit mounting
Features #3	Cable entry slots underside of housing