



## RTD Miniatur 3-poliger Anschlussstecker



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | XE-7770-001                            |
| <b>General Description</b> | Miniature connectors for use with RTDs |
| <b>Connector Size</b>      | Miniature                              |
| <b>Connector Type</b>      | In-Line Plug                           |
| <b>Material (Body)</b>     | Glass Filled Thermoplastic             |
| <b>Material (Pins)</b>     | Copper (Cu)                            |
| <b>Colour</b>              | White                                  |
| <b>Colour Code</b>         | IEC 60584.3                            |
| <b>Max. Temperature</b>    | 220°C                                  |
| <b>Pins</b>                | Flat 3-Pins                            |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                 |