


RTD-Miniatur-3-polige Steckverbinderbuchse



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Product Code	XE-7771-001
General Description	Miniature connectors for use with RTDs
Connector Size	Miniature
Connector Type	In-Line Socket
Material (Body)	Glass Filled Thermoplastic
Material (Pins)	Copper (Cu)
Colour	White
Colour Code	IEC 60584.3
Max. Temperature	220°C
Pins	Accepts Flat 3-Pins
Standards Met	BS EN 50212:1997. RoHS