


Miniatur-Schnelldraht-Thermoelement-Steckverbinderbuchse AM-K-FQ Typ K ANSI



Miniature in-line Socket - Quick Wire
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Schnellverdrahtungs-Thermoelement-Buchse AM-K-FQ Typ K ANSI

Einfache Jab-In-Verbindung, einfach den Draht eindrcken und die Schraube festziehen. Die Kontakte sind polarisiert, um einen falschen Anschluss zu verhindern.

Leistungsbeschreibung**Specifications**

Product Code	XE-1142-001
General Description	Easy jab-in connection, just push in wire & tighten screw.
Thermocouple Type	K ANSI
Connector Type	Quick Wire Socket
Connector Size	Miniature
Colour	Yellow
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector