

WT9055 Digitaler Schallpegelmesser



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Der digitale Schallpegelmesser WT9055 ist neuartig, klein und tragbar. Das Gert eignet sich fr die Messung von Fabrikumgebungen, Qualittskontrolle, Gesundheitsprvention und verschiedene Umgebungsgerusche, einschlielich Bros, Transportwege, Stereoanlagen und vieles mehr.

Messung des Schallpegels (dBA)

MIN / MAX / Stromwert sperren

Daten-Hold

LCD-Hintergrundbeleuchtungsfunktion

Manuelle / automatische Abschaltung

Alarmierung bei Hintergrundbeleuchtung

Leistungsbeschreibung

Specifications

Product Code	XE-4157-001
Model	WT9055
Type	Sound Level Meter
Accuracy	+/- 1.5dB
Resolution	0.1dB(A)
Battery	3 x 1.5V AAA
Dimensions	53.4mm x 31.4mm x 131.7mm
Working Temperature	0-40°C
Operating Humidity	10-80% RH
Weight	96.38g
Storage Temperature	10-60°C
Measuring range	35-130dB (A)
Important Information	Frequency response: 31.5Hz - 8KHz