


# Miniatur-Schnelldraht-Thermoelement-Steckverbinderbuchse AM-N-FQ Typ N ANSI



**Miniature in-line Socket - Quick Wire**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur-Schnellkupplungs-Thermoelement-Buchse AM-N-FQ Typ N ANSI**

Einfache Jab-In-Verbindung, einfach den Draht eindrcken und die Schraube festziehen. Die Kontakte sind polarisiert, um einen falschen Anschluss zu verhindern.

**Leistungsbeschreibung****Specifications**

<b>Product Code</b>	XE-1282-001
<b>General Description</b>	Easy jab-in connection, just push in wire & tighten screw.
<b>Thermocouple Type</b>	N ANSI
<b>Connector Type</b>	Quick Wire Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Orange
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel / Chromium / Silicon
<b>-Negative Contact</b>	Nickel Silicon
<b>Standards Met</b>	BS EN 50212:1997. RoHS