

# Miniatur-Thermoelement-Steckverbinder mit rechteckiger Blende AM-K-FF Typ K ANSI



**Miniature Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur Rechteck Thermoelement Steckverbinder Blende Buchse AM-K-FF Typ K ANSI**

Nylonclip für eine schnelle und einfache Befestigung an der Platte. Alle Blenden-Montagesockel sind mittels ungleich großer Stifte polarisiert

**Leistungsbeschreibung****Specifications**

<b>Product Code</b>	XE-1137-001
<b>General Description</b>	All fascia mounting sockets are polarised by means of unequally sized pins
<b>Thermocouple Type</b>	K ANSI
<b>Connector Type</b>	Fascia mounting socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Yellow
<b>Dimensions</b>	See Schematic Drawing (Product Literature Tab)
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS
<b>Series</b>	Miniature Connector