


Miniatur-Thermoelement-Steckverbinder mit rechteckiger Blendenbuchse AM-J-FF Typ J ANSI



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Rechteckiger Thermoelement-Steckverbinder Blendenbuchse AM-J-FF Typ J ANSI

Nylon-Clip zur schnellen und einfachen Befestigung an der Platte. Alle Blendenmontagebuchsen sind mittels ungleich groer Stifte polarisiert

Leistungsbeschreibung**Specifications**

| | |
|----------------------------|--|
| Product Code | XE-9067-001 |
| General Description | All fascia mounting sockets are polarised by means of unequally sized pins |
| Thermocouple Type | J ANSI |
| Connector Type | Fascia mounting socket |
| Connector Size | Miniature |
| Colour | Black |
| Max. Temperature | 120°C |
| +Positive Contact | Iron |
| -Negative Contact | Copper Nickel |
| Standards Met | BS EN 50212:1997. RoHS |