

# Miniatur Rechteckige Fascia Socket Thermoelement Stecker AM-R/S-FF Typ R/S ANSI



**Miniature Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur Rechteckige Thermoelement Stecker Fascia Socket AM-R/S-FF Typ R/S ANSI**

Nylon-Clip für schnelle, einfache Befestigung am Panel. Alle Faszien-Montagebuchsen werden durch ungleich große Stifte polarisiert

**Spezifikationen****Specifications**

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | XE-1317-001  |
| <b>General Description</b> | All fascia mounting sockets are polarised by means of unequally sized pins |
| <b>Thermocouple Type</b>   | R/S ANSI   |
| <b>Connector Type</b>      | Fascia mounting socket   |
| <b>Connector Size</b>      | Miniature  |
| <b>Colour</b>              | Green  |
| <b>Max. Temperature</b>    | 120°C  |
| <b>+Positive Contact</b>   | Copper   |
| <b>-Negative Contact</b>   | Copper Alloy   |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS   |