

Miniatur-Panel-Montage Thermoelement-Steckverbinder mit Edelsthalhalterung AM-E-SSPF Typ E ANSI



Miniature in-line Socket with stainless steel mounting bracket
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Panel-Montage Thermoelement-Steckverbinder mit Edelstahlhalterung AM-E-SSPF Typ E ANSI

Edelstahlhalterung einfach am Stecker montiert. Vollständig geschlossene Klemmen mit Kabelanschlussdichtung.

Spezifikationen**Specifications**

Product Code	XE-9038-001
General Description	Fully enclosed terminals with cable port seal
Thermocouple Type	E ANSI
Connector Type	Socket with SS Bracket
Connector Size	Miniature
Colour	Violet
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS