




Leiterplattenmontage Thermoelement-Steckverbinderbuchse AM-K-PCB Typ K ANSI



Miniature PCB Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leiterplattenmontage Thermoelement-Steckverbinderbuchse AM-K-PCB Typ K ANSI

Speziell fr Leiterplatten, massive Pins fr die direkte Platinenmontage und Bereitstellung fr CJC-Sensor

Leistungsbeschreibung**Specifications**

| | |
|----------------------------|---|
| Product Code | XE-1139-001 |
| General Description | Versatile, high quality circuit board mounting miniature socket |
| Thermocouple Type | K ANSI |
| Connector Type | PCB socket |
| Connector Size | Miniature |
| Colour | Yellow |
| Max. Temperature | 220°C |
| +Positive Contact | Nickel Chromium |
| -Negative Contact | Nickel Alloy |
| Standards Met | BS EN 50212:1997. RoHS |