


Standard-Thermoelement-Anschlussstecker AS-K-M Typ K ANSI



Standard in-line Plug
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Thermoelement-Anschlussstecker AS-K-M Typ K ANSI

Die Kontakte bestehen aus echten Thermoelement-Legierungen. Geeignet für industrielle Anwendungen

Leistungsbeschreibung**Specifications**

Product Code	XE-1430-001
General Description	Contacts are polarised to prevent incorrect connection
Thermocouple Type	K ANSI
Connector Type	Plug
Connector Size	Standard
Colour	Yellow
Dimensions	See Schematic Drawing (Product Literature Tab)
Max. Temperature	220°C
Pins	Solid Round Pins
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS