



Standard Thermoelement-Steckverbinder In-Line-Sockel AS-N-F Typ N ANSI



Standard in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard Thermoelement-Steckverbinder In-Line-Sockel AS-N-F Typ N ANSI

Akzeptiert Standard-Rundstecker, geeignet für industrielle Anwendungen

Spezifikationen**Specifications**

Product Code	XE-1572-001
General Description	Accepts standard round plug pins, contacts are polarised to prevent incorrect connection.
Thermocouple Type	N ANSI
Connector Type	Socket
Connector Size	Standard
Colour	Orange
Max. Temperature	220°C
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS