


# Rechteckige Standard-Blenden-Thermoelement-Steckverbinde rbuchse AS-K-FF Typ K ANSI



**Standard Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Standard-Rechteck-Thermoelement-Steckverbinder Blendenbuchse AS-K-FF Typ K ANSI**

Nylon-Clip für eine schnelle und einfache Befestigung an der Platte. Alle Fasziensbefestigungsbuchsen sind mittels ungleich großer Stifte polarisiert

**Leistungsbeschreibung****Specifications**

<b>Product Code</b>	XE-1437-001
<b>General Description</b>	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
<b>Thermocouple Type</b>	K ANSI
<b>Connector Type</b>	Fascia Mounting Socket
<b>Connector Size</b>	Standard
<b>Colour</b>	Yellow
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS