



Standard-rechteckige Blenden-Thermoelement-Steckverbinderbuchse AS-T-FF Typ T ANSI



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Rechteckiger Thermoelement-Steckverbinder Blendenbuchse AS-T-FF Typ T ANSI

Nylon-Clip zur schnellen und einfachen Befestigung an der Platte. Alle Blendenmontagebuchsen sind mittels ungleich groer Stifte polarisiert

Leistungsbeschreibung**Specifications**

Product Code	XE-1527-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	T ANSI
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Blue
Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS