



# Standard-rechteckige Blenden-Thermoelement-Steckverbinderbuchse AS-T-FF Typ T ANSI



**Standard Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Standard-Rechteckiger Thermoelement-Steckverbinder Blendenbuchse AS-T-FF Typ T ANSI**

Nylon-Clip zur schnellen und einfachen Befestigung an der Platte. Alle Blendenmontagebuchsen sind mittels ungleich groer Stifte polarisiert

**Leistungsbeschreibung****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-1527-001   |
| <b>General Description</b> | Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel |
| <b>Thermocouple Type</b>   | T ANSI  |
| <b>Connector Type</b>      | Fascia Mounting Socket  |
| <b>Connector Size</b>      | Standard  |
| <b>Colour</b>              | Blue  |
| <b>Max. Temperature</b>    | 120°C   |
| <b>+Positive Contact</b>   | Copper  |
| <b>-Negative Contact</b>   | Copper Nickel   |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS  |