




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Standard Rechteckige Fascia Thermoelement Stecker Socket AS-N-FF Typ N ANSI



Standard Fascia Socket (with nylon mounting clip)

(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard Rechteckiger Thermoelement-Steckverbinder Fascia Socket AS-N-FF Typ N ANSI

Nylon-Clip für schnelle, einfache Befestigung am Panel. Alle Faszien-Montagebuchsen werden durch ungleich große Stifte polarisiert

Spezifikationen**Specifications**

Product Code	XE-1577-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	N ANSI
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Orange
Max. Temperature	120°C
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS