


# Thermoelement-Steckverbinder für die Standard-Panelhalterung mit Edelstahlhalterung AS-R/S-SSPF Typ R/S ANSI



**Standard in-line Socket with stainless steel mounting bracket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Thermoelement-Steckverbinder für die Standard-Panelhalterung mit Edelstahlhalterung AS-R/S-SSPF Typ  
R/S ANSI**

Für den Einsatz in industriellen Prozessen vorgesehen, kann aber in allen Anwendungen nach Belieben eingesetzt werden

**Spezifikationen****Specifications**

<b>Product Code</b>	XE-1618-001
<b>General Description</b>	Intended for use in more industrial processes but can be used in all applications as desired
<b>Thermocouple Type</b>	R/S ANSI
<b>Connector Type</b>	Socket with SS Bracket
<b>Connector Size</b>	Standard
<b>Colour</b>	Green
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS