



# Miniatur-Schnelldraht Thermoelement Stecker JM-K-MQ Typ K JIS



**Miniature in-line Plug - Quick Wire**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur Thermoelement Stecker SchnelldrahtStecker JM-K-MQ Typ K JIS**

Miniaturstecker mit Flächensteckern. Einfacher Jab-in-Anschluss für schnelle Beendigung. Einfach einschieben Draht & Ziehen Sie Schraube.

**Spezifikationen****Specifications**

<b>Product Code</b>	XE-9162-001
<b>General Description</b>	Easy jab-in connection for rapid termination.
<b>Thermocouple Type</b>	K JIS
<b>Connector Type</b>	Quick Wire Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	Blue
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Flat Pins
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS