


Miniatur-Schnelldraht-Thermoelement-Steckverbinderbuchse JM-J-FQ Typ J JIS



Miniature in-line Socket - Quick Wire
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Schnellkupplungs-Thermoelement-Buchse JM-J-FQ Typ J JIS

Einfache Jab-In-Verbindung, einfach den Draht eindrcken und die Schraube festziehen. Die Kontakte sind polarisiert, um einen falschen Anschluss zu verhindern.

Leistungsbeschreibung**Specifications**

Product Code	XE-9183-001
General Description	Easy jab-in connection, just push in wire & tighten screw.
Thermocouple Type	J JIS
Connector Type	Quick Wire Socket
Connector Size	Miniature
Colour	Yellow
Max. Temperature	220°C
+Positive Contact	Iron
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS