

Miniatur Rechteckige Faszien Sockel Thermoelement Stecker JM-E-FF Typ E JIS



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur Rechteckige Thermoelement Stecker Fascia Sockel JM-E-FF Typ E JIS

Nylon-Clip für schnelle, einfache Befestigung am Panel. Alle Faszien-Montagebuchsen werden durch ungleich große Stifte polarisiert

Spezifikationen**Specifications**

| | |
|----------------------------|----------------------------------------------------------------------------|
| Product Code | XE-9237-001 |
| General Description | All fascia mounting sockets are polarised by means of unequally sized pins |
| Thermocouple Type | E JIS |
| Connector Type | Fascia mounting socket |
| Connector Size | Miniature |
| Colour | Violet |
| Max. Temperature | 120°C |
| +Positive Contact | Nickel Chromium |
| -Negative Contact | Copper Nickel |
| Standards Met | BS EN 50212:1997. RoHS |