

Miniatur-Panel-Montage Thermoelement-Steckverbinder mit Edelsthalterung JM-J-SSPF Typ J JIS



Miniature in-line Socket with stainless steel mounting bracket
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Panel-Montage Thermoelement-Steckverbinder mit Edelstahlhalterung JM-J-SSPF Typ J JIS

Edelstahlhalterung einfach am Stecker montiert. Vollständig geschlossene Klemmen mit Kabelanschlussdichtung

Spezifikationen**Specifications**

Product Code	XE-9178-001
General Description	Fully enclosed terminals with cable port seal
Thermocouple Type	J JIS
Connector Type	Socket with SS Bracket
Connector Size	Miniature
Colour	Yellow
Max. Temperature	220°C
+Positive Contact	Iron
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS