


# Miniatur-Panel-Montage Thermoelement-Steckverbinder mit Edelstahlhalterung JM-B-SSPF Typ B JIS



**Miniature in-line Socket with stainless steel mounting bracket**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur-Panel-Montage Thermoelement-Steckverbinder mit Edelstahlhalterung JM-B-SSPF Typ B JIS**

Edelstahlhalterung einfach am Stecker montiert. Vollständig geschlossene Klemmen mit Kabelanschlussdichtung

**Spezifikationen****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-9218-001                                   |
| <b>General Description</b> | Fully enclosed terminals with cable port seal |
| <b>Thermocouple Type</b>   | B JIS   |
| <b>Connector Type</b>      | Socket with SS Bracket                        |
| <b>Connector Size</b>      | Miniature                                     |
| <b>Colour</b>              | Grey  |
| <b>Max. Temperature</b>    | 220°C   |
| <b>+Positive Contact</b>   | Copper  |
| <b>-Negative Contact</b>   | Copper  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                        |