


Standard Thermoelement Connector In-Line Socket JS-K-F Typ K JIS



Standard in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard Thermoelement Connector In-Line Socket JS-K-F Typ K JIS

Akzeptiert Standard-Rundstecker, geeignet für industrielle Anwendungen

Spezifikationen**Specifications**

Product Code	XE-9272-001
General Description	Accepts standard round plug pins, contacts are polarised to prevent incorrect connection.
Thermocouple Type	K JIS
Connector Type	Socket
Connector Size	Standard
Colour	Blue
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS