





# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Standard Rechteckige Fascia Thermoelement Stecker Sockel JS-R/S-FF Typ R/S JIS



**Standard Fascia Socket (with nylon mounting clip)**

(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Standard Rechteckiger Thermoelement-Steckverbinder Fascia Socket JS-R/S-FF Typ R/S JIS**

Nylon-Clip für schnelle, einfache Befestigung am Panel. Alle Faszien-Montagebuchsen werden durch ungleich große Stifte polarisiert

**Spezifikationen****Specifications**

<b>Product Code</b>	XE-9317-001
<b>General Description</b>	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
<b>Thermocouple Type</b>	R/S JIS
<b>Connector Type</b>	Fascia Mounting Socket
<b>Connector Size</b>	Standard
<b>Colour</b>	Black
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS