





LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Standard Rechteckige Fascia Thermoelement Stecker Sockel JS-R/S-FF Typ R/S JIS



Standard Fascia Socket (with nylon mounting clip)

(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard Rechteckiger Thermoelement-Steckverbinder Fascia Socket JS-R/S-FF Typ R/S JIS

Nylon-Clip für schnelle, einfache Befestigung am Panel. Alle Faszien-Montagebuchsen werden durch ungleich große Stifte polarisiert

Spezifikationen**Specifications**

Product Code	XE-9317-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	R/S JIS
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Black
Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper Alloy
Standards Met	BS EN 50212:1997. RoHS