

Barriereklemme ThermoelementStecker JS-J-BTS Typ J JIS



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Barriereklemme ThermoelementStecker JS-J-BTS Typ J JIS

Lieferung in einzelnen, ineinandergreifenden Paaren. Anschlüsse im Thermoelementmaterial

Spezifikationen**Specifications**

Product Code	XE-1961-001
General Description	Supplied in single, interlocking pairs. Connections in thermocouple material
Thermocouple Type	J JIS
Connector Type	Barrier Terminal Blocks
Colour	Yellow
Dimensions	See Schematic Drawing (Product Literature Tab)
Max. Temperature	220°C
+Positive Contact	Iron
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS