



**LABFACILITY**

TEMPERATURE & PROCESS TECHNOLOGY

## Typ K Einziehbares curly Thermoelement Lead (IEC)




XE-7030-001

**TYPE-K (IEC)**



XE-7030-001

**TYPE-K (IEC)**



XE-7030-001

**TYPE-K (IEC)**

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Verfügbar in Typ K

Farbe kodiert nach IEC

2 Meter lang bei vollständiger Verlängerung

Beendet mit einem eingebauten Miniaturstecker

Polyurethan außen

## Specifications

<b>Product Code</b>	XE-7030-001
<b>Cable Length</b>	2M fully extended
<b>Cable Type</b>	Polyurethane outer
<b>Connector Type</b>	Thermocouple Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	Green
<b>Standards Met</b>	Colour code (cores & jacket) to IEC-584-3
<b>Thermocouple Type</b>	K
<b>Termination Type</b>	Miniature Thermocouple Plug