




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Typ K Einziehbares curly Thermoelement Lead (IEC)




XE-7030-001

TYPE-K (IEC)



XE-7030-001

TYPE-K (IEC)



XE-7030-001

TYPE-K (IEC)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Verfügbar in Typ K

Farbe kodiert nach IEC

2 Meter lang bei vollständiger Verlängerung

Beendet mit einem eingebauten Miniaturstecker

Polyurethan außen

Specifications

Product Code	XE-7030-001
Cable Length	2M fully extended
Cable Type	Polyurethane outer
Connector Type	Thermocouple Plug
Connector Size	Miniature
Colour	Green
Standards Met	Colour code (cores & jacket) to IEC-584-3
Thermocouple Type	K
Termination Type	Miniature Thermocouple Plug